

# Synthesis of Metal Colloids

---

**Caution:** Since the colloid synthesis is sensitive to various laboratory conditions and the practices of the individual researchers, the following procedure should be taken as a general guideline. One may have to vary the experimental conditions to suit their needs.

## Preparation of $2.5 \times 10^{-4}$ M Gold Colloids (Sodium Citrate Reduction Method)

1. Make a solution of  $\sim 5.0 \times 10^{-3}$  M  $\text{HAuCl}_4$  in water. (0.1699 g  $\text{HAuCl}_4$  in 100 mL deionized  $\text{H}_2\text{O}$ )
2. Take 1 mL of that solution and add it to another 18 mL of  $\text{H}_2\text{O}$ .
3. Make a solution of 0.5% sodium citrate (0.25g in 50 mL of  $\text{H}_2\text{O}$ ).
4. Heat the 19 mL solution of  $\text{HAuCl}_4$  until it begins to boil.
5. Add 1 mL of 0.5% sodium citrate solution, as soon as boiling commences.
6. Continue heating until colour change is evident (pale purple).
7. Remove the solution from the heating element and continue to stir until it has cooled to room temperature.
8. Top the solution up to 20 mL to account for boiling.

*This is an old method (referred as Turkevich method) which yields fairly uniform size colloids with diameter of 15-20 nm. See reference Turkevich, J.;Stevenson, P. L.;Hillier, J.Discuss. Faraday Soc. 1951, 11, 55*

---

## Preparation of Gold Particles in Toluene:

- A. Hydrogen tetrachloroaurate (30 mL of 30 mM, in water).

$$\begin{aligned}\text{Mass} &= 393\text{g/mol} \times 0.03 \text{ M} \times 30\text{mL}/1000\text{mL} \\ &= 0.3537\text{g of HAuCl}_4 \text{ in 30mL of H}_2\text{O}\end{aligned}$$

- B. Tetraoctyl ammonium bromide (80 mL in 50 mM, in toluene).

$$\begin{aligned}\text{Mass} &= 546.8\text{g/mol} \times 0.05 \text{ M} \times 80\text{ml}/1000\text{mL} \\ &= 2.187\text{g of TOAB in 80mL of toluene}\end{aligned}$$

1. Prepare 2.19 g of tetraoctyl ammonium bromide in 80 mL of toluene.

2. Add solution prepared in step 1. to a solution of hydrogen tetrachloroaurate (0.3537 g in 30 mL of H<sub>2</sub>O).
3. Stir for 10 min.
4. Vigorously stir reaction mixture and add NaBH<sub>4</sub> (0.38 g in 25 mL of H<sub>2</sub>O) dropwise over a period of ~30 min. (Ensure that organic and aqueous phases are being mixed together).
5. Stir solution for an additional 20 min.
6. Extract organic phase and wash once with diluted H<sub>2</sub>SO<sub>4</sub> (for neutralization) and five times with distilled water.
7. Dry organic layer with Na<sub>2</sub>SO<sub>4</sub>.

For platinum particles substitute dihydrogen hexachloroplatinate (IV) (PtCl<sub>6</sub> H<sub>2</sub>O) for HAuCl<sub>4</sub> and for iridium particles substitute H<sub>2</sub>IrCl<sub>6</sub>·4H<sub>2</sub>O for HAuCl<sub>4</sub>.

Yields highly concentrated gold colloidal suspension with particle diameter in the range of 5-10 nm. Can be suspended in both polar and nonpolar solvents.

Adopted from the reference Brust, M.; Walker, M.; Bethell, D.; Schiffrin, D. J. Whyman, R., *J. Chem. Soc., Chem. Commun.*, 1994, 801-802 and George Thomas, K. Kamat, P. V. *J. Am. Chem. Soc.* 2000, 122, 2655

---

## Preparation of 1.0 x 10<sup>-3</sup>M Ag Colloids (Sodium Citrate Reduction Method)

1. Make a solution of ~ 5.0 x 10<sup>-3</sup>M AgNO<sub>3</sub> in water. (0.0425 g in 50 mL deionized H<sub>2</sub>O).
2. Take 25 mL of that solution and add it to another 100 mL of H<sub>2</sub>O (now ~1.0 x 10<sup>-3</sup>M).
3. Make a solution of 1% sodium citrate (0.5 g in 50 mL of H<sub>2</sub>O).
4. Heat the 125 mL solution of AgNO<sub>3</sub> until it begins to boil.
5. Add 5 mL of 1% sodium citrate solution, as soon as boiling commences.
6. Continue heating until a colour change is evident (pale yellow).
7. Remove the solution from the heating element and continue to stir until it has cooled to room temperature.
8. Top the solution up to 125 mL to account for boiling.

This method yields relatively large size silver nanocrystallites with a diameter of 60-80 nm and exhibits abs. max. ~420 nm. See reference:

*J. Phys. Chem. B*, 1998, 102, 3123

*(Note: Use of Sodium Borohydride as a reductant can give smaller size silver nanoparticles with plasmon absorption around 380 nm. Presence of citric acid or polyvinyl alcohol can provide additional stability to these colloids)*

---

## Preparation of Au Capped Ag

1. Take 125 mL of  $1.0 \times 10^{-3}$  M solution of Ag colloids and add 12,5 mL H<sub>2</sub>O.
  2. Heat this solution until it comes to a boil and then add the appropriate amount of  $5.0 \times 10^{-3}$  M HAuCl<sub>4</sub>. (For example: 500 μL, 1000 μL, 1500 μL, 3000 μL, and 5000 μL)
- 

## Silver Nanocrystal Preparation in Organic Medium (Using TOAB)

Modification of method found in source: Korgel, B.A.; Fullam, S.; Connolly, S.; Fitzmaurice, D. *J. Phys. Chem. B* **1998**, *102*, 8379-8388. (The method described in this paper yields a precipitate of AgBr in the initial extraction process)

1. Prepare ~5.0M NaNO<sub>3</sub> in deionized water (12.749g NaNO<sub>3</sub> in 30mL H<sub>2</sub>O).
2. Prepare ~50mM TOAB in toluene (1.367g tetraoctylammonium bromide in 50mL toluene).
3. Add the TOAB/toluene solution to the NaNO<sub>3</sub>/water solution.
4. Stir vigorously for 1 hour (to remove Br<sup>-</sup> ions from solution and prevent the formation of AgBr when AgNO<sub>3</sub> is added).
5. Extract organic phase and set aside. Discard aqueous phase.
6. Prepare ~30mM AgNO<sub>3</sub> in water (0.0764g NaNO<sub>3</sub> in 14mL H<sub>2</sub>O).
7. Add 7.5mL of 30mM AgNO<sub>3</sub> solution to the organic solution.
8. Stir vigorously for 45 minutes.
9. Extract organic phase (discard aqueous layer).
10. Add 0.16mg (~0.189mL) of 1-dodecanethiol to organic solution (to cap the silver)
11. Stir vigorously for 15 minutes.
12. Meanwhile, prepare ~0.4M NaBH<sub>4</sub> in water (0.3783g NaBH<sub>4</sub> in 24mL H<sub>2</sub>O).
13. Add 6.25mL of the NaBH<sub>4</sub>, dropwise over a 35min. period, to the solution containing the silver (organic layer), while stirring vigorously.
14. Stir for ~15 hours (overnight).
15. Extract organic layer (discard aqueous layer).
16. Wash organic layer 3 times with dilute ethanol.
17. Allow to settle, and extract organic layer.
18. Store in closed container.